



# Information Sheet for BioMedicure's Products

---

## Vaccination Schedule:

Day 1 and day 15

The first vaccination is mostly for antigen presentation and memory and polyclonal antibody production. The second vaccination will synergize immune components including polyclonal antibody, lymphocytes, interleukin and various cytokines to mount attacks against multiple types of cancerous cells including pre-cancer cells, cancer cells, metastasized cancer cells and cancer stem cells.

## Vaccination Protocol

1. Read user's manual for instructions.
2. Label unvaccinated mice with ear punch or other methods.
3. Thaw the vial by a gentle agitation. Thawing should be rapid (within half an hour).
4. Using 1-2 mL syringes with clear marks on 100 uL and a 25-27 gauge needle.
5. The site of vaccination may be shaved so the sub-Q vaccination is not interfered with by fur.
6. Inject 100 uL CancerVaccine™ by sub-Q.
7. Push syringe slowly (30 seconds for 100 uL) to give optimal sub-Q presentation.

## Additional Information

Additional product and technical information can be obtained from the catalog references and the BioMedicure Web site at [www.biomedicure.com](http://www.biomedicure.com), or by email at [cv@biomedicure.com](mailto:cv@biomedicure.com)

## BioMedicure Warranty

The efficacy of BioMedicure products is warranted for 30 days from the date of shipment. If you feel there is a problem with this product, contact Technical Services by phone at 858-586-1888 or by email at [cv@biomedicure.com](mailto:cv@biomedicure.com).

## Disclaimers

This product is intended for laboratory research purposes only. It is not intended for use in any regulatory area or subjects including human beings.

While BioMedicure uses reasonable efforts to include accurate and up to date information on this product sheet, BioMedicure makes no warranties or representations as to its accuracy. Citations from scientific literature and patents are provided for informational purposes only. BioMedicure does not warrant that such information has been confirmed to be accurate.

This product is sent with the condition that you are responsible for its safe storage, handling and use. BioMedicure is not liable for any damages or injuries arising for receipt and/or use of this product. While reasonable effort is made to insure authenticity and reliability of strains on deposit, BioMedicure is not liable for damages arising from the misidentification or misrepresentation of cultures.

---

**BioMedicure**

7940 Silverton Avenue, Suite 107  
San Diego, CA 92126

Tel: 858-586-1888  
Fax: 858-997-2598  
Email: [info@biomedicure.com](mailto:info@biomedicure.com)  
[www.biomedicure.com](http://www.biomedicure.com)