

Information Sheet for BioMedicure's Products

Rat Urinary Bladder Tumor CancerVaccine™

BioMedicure Number:	CV-655-5	Price:	\$1299.99
Designations:	V-655		
Unit Size:	1000 µL		
Producer:	BioMedicure		
Biosafety level:	0		
Shipped:	Frozen	Use Restrictions	These cells are distributed for research purposes only. BioMedicure recommends that individuals contemplating commercial use of any enriched CancerStemCell™ products first contact BioMedicure for potential licensing agreement. Buyers are the second party end users of this product. Any distribution of this product to a third party is strictly prohibited.
Organism:	<i>Rattus norvegicus</i>		
Strain:	Wistar		
Target disease:	Tumor, urinary bladder		
Growth Properties:	Non-viable	Biosafety Level: 0	Appropriate safety procedures should always be used with this material. Laboratory safety is discussed in the following publications: <i>Biosafety in Microbiological and Biomedical Laboratories</i> , 4 th ed. HHS Publication No. (CDC) 93-8395. U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Washington DC: U.S. Government Printing Office: 1999. The entire text is available online at www.cdc.gov/od/ohs/biosfty/bmbl4/bmbl4toc.htm
Tumorigenic <i>in vitro</i>:	No		
Tumorigenic <i>in vivo</i>:	No		
Toxicity:	None		
Heterogeneous:	Yes		
Preventive Cancer:	Tumor, urinary bladder	Handling Procedure for Frozen CancerVaccine™	To insure the highest level of efficacy, thaw the vial and initiate the vaccination within half an hour. If upon arrival, continued storage of the frozen CancerVaccine™ is necessary, it should be stored at -20°C. Storage at -70°C or lower temperature is not required. Storage at room temperature or 4-8°C not recommended as will result in antigen loss leading to reduced efficacy in immune response.
Paired CancerStemCell™:	CSC-655		
Molecular recognition patterns:	None		
Genomic background:	Rat Wistar		
Antigens:	Diverse		
Immune response rate:	100%	SAFETY PRECAUTION:	BioMedicure highly recommends that protective gloves and clothing always be used.

BioMedicure

7940 Silvertown Avenue, Suite 107
San Diego, CA 92126

Tel: 858-586-1888
Fax: 858-997-2598
Email: info@biomedicure.com
www.biomedicure.com

Information Sheet for BioMedicure's Products

Vaccination Schedule:

Day 1 and day 15

The first vaccination is mostly for antigen presentation and memory and polyclonal antibody production. The second vaccination will synergize immune components including polyclonal antibody, lymphocytes, interleukin and various cytokines to mount attacks against multiple types of cancerous cells including pre-cancer cells, cancer cells, metastasized cancer cells and cancer stem cells.

Vaccination Protocol

1. Read user's manual for instructions.
2. Label unvaccinated mice with ear punch or other methods.
3. Thaw the vial by a gentle agitation. Thawing should be rapid (within half an hour).
4. Using 1-2 mL syringes with clear marks on 100 uL and a 25-27 gauge needle.
5. The site of vaccination may be shaved so the sub-Q vaccination is not interfered with by fur.
6. Inject 100 uL CancerVaccine™ by sub-Q.
7. Push syringe slowly (30 seconds for 100 uL) to give optimal sub-Q presentation.

Additional Information

Additional product and technical information can be obtained from the catalog references and the BioMedicure Web site at www.biomedicure.com, or by email at cv@biomedicure.com

BioMedicure Warranty

The efficacy of BioMedicure products is warranted for 30 days from the date of shipment. If you feel there is a problem with this product, contact Technical Services by phone at 858-586-1888 or by email at cv@biomedicure.com.

Disclaimers

This product is intended for laboratory research purposes only. It is not intended for use in any regulatory area or subjects including human beings.

While BioMedicure uses reasonable efforts to include accurate and up to date information on this product sheet, BioMedicure makes no warranties or representations as to its accuracy. Citations from scientific literature and patents are provided for informational purposes only. BioMedicure does not warrant that such information has been confirmed to be accurate.

This product is sent with the condition that you are responsible for its safe storage, handling and use. BioMedicure is not liable for any damages or injuries arising for receipt and/or use of this product. While reasonable effort is made to insure authenticity and reliability of strains on deposit, BioMedicure is not liable for damages arising from the misidentification or misrepresentation of cultures.

BioMedicure

7940 Silverton Avenue, Suite 107
San Diego, CA 92126

Tel: 858-586-1888
Fax: 858-997-2598
Email: info@biomedicure.com
www.biomedicure.com